

Technical Data Sheet

SafeBasics[™] Saliva Ejectors - Clear Blue Tip Item# SFBGSE0782

Item Code & Color Item Code Color

SFBGSE0782 Clear cannula with blue fixed tip

Intended Use Saliva ejectors are used in dentistry in order to remove saliva,

blood or other debris from the mouth during dental procedures.

This disposable device is

Packaging 100 units per bag

10 bags per case 1000 units per case

Material PVC Cannula, reinforced with zinc or copper

covered wire or anneaded steel wire core.


Size Length of Saliva Ejector 145 -150 mm

Diameter of Cannula 6.5 – 6.6 mm

Usage For single-use only

Sterility Non-Sterile

Country of Origin Italy



This technical data sheet issued by A.R. Medicom Inc (Asia) Ltd

TDS-SFBGSE0782 Revision 00 Effective Date: 2015-07-06 Page 1 of 1